


IntelConnect Setup Guide

STEP 1 – Set up IntelConnect

1. Open your web browser
2. In the address bar, enter your PACS link:


All other Canopy Imaging Regions

<https://trgpacs.co.nz/Intelconnect>



Manawatu & Horowhenua Region

<https://images.pnxray.co.nz/IntelConnect>



3. You will land on the IntelConnect login page
4. If this is your first time logging-in, you'll be prompted to activate Multi-Factor Authentication (MFA); your MFA set-up guide attached to this email
5. Log in with your Canopy Imaging PACS username and password
6. Once you accept the confidentiality agreement, you've completed your MFA and IntelConnect setup.

STEP 2 – Navigating IntelConnect to view images and reports


Click on the Profile tab and review your details, if anything needs updating, please email us:

referrersupport@canopyimaging.co.nz



STEP 3 – Searching for imaging and reports (if a break glass button appears next to the search function see step 3A at the bottom of this guide)

- To access imaging or reports for your patients or patients coming into your care (that were not originally referred by you):
 - Start typing the NHI or name in the tab 'Search patient name, id, accession'
 - Click the magnifying glass icon to search



If you see this warning about your RIS ID, you can ignore and still search for patient imaging and reports

 Your user account does not have a RIS ID. Referred patients cannot be found without a RIS ID. Please contact your radiology provider and ask for a RIS ID to be added to your user account.

Note: a RIS ID is required to be mapped to your IntelConnect PACS account to setup notifications, if you'd like us to set this up, please fill out our new referrer form here: [new referrer form](#) and we'll email once this has been actioned

STEP 4 – View Images and Reports

- When images and/or reports are available, you'll see the image icon and report icon in the search results

kestral Search Refine your search with Advanced Search						
Kestral, Hamish (test) DOB: 01-Jan-2001 (25 Y) Gender: Male Patient ID: PUA0569 125 orders Request an exam						
Accession Number	Mod.	Exam Description	Exam Date	Exam Details	Org	Image Availability
2026AH0003504	MR	MRI Right Shoulder	12-Aug-2026 10:51	 Order Complete (CM)	AH	NO IMAGES
2026NP0004031	CR	XR Right Shoulder	12-Aug-2026 10:44	 Order Complete (CM)	NP	NO IMAGES
2026AH0003503	US	US Right Shoulder	12-Aug-2026 10:42	 Order Complete (CM)	AH	NO IMAGES
2026NP0004020	CR	XR Right Knee	11-Aug-2026 21:40	 Final Report (ZZ)	NP	NO IMAGES
2026AH0003053	CR	XR Lumbar Spine	13-Jul-2026 14:36	 Final Report (ZZ)	AH	NO IMAGES
2026AH0003052	CR	XR Left Shoulder	13-Jul-2026 14:31	 Final Report (ZZ)	AH	NO IMAGES
2026HS0007055	CR	XR Lumbar Spine	29-Jun-2026 15:46	 Final Report (ZZ)	HS	NO IMAGES
2026NP0003283	US	US Abdomen	29-Jun-2026 15:40	 Final Report (ZZ)	NP	NO IMAGES
2026NP0003274	CR	XR Right Shoulder	29-Jun-2026 14:04	 Final Report (ZZ)	NP	NO IMAGES
2026LI0007594	CR	XR Left Finger Hand and Wrist	19-May-2026 23:12	 Final Report (ZZ)	LI	NO IMAGES
2026LI0006466	US	US Abdomen and Pelvis	01-May-2026 15:07	 Final Report (ZZ)	LI	ONLINE  
2026SH0004044	CR	XR Left Ankle	12-Mar-2026 12:40	 Final Report (ZZ)	SH	ONLINE  
2025NX0015035	US	US Abdomen	17-Nov-2025 17:28	 Final Report (ZZ)	NX	ONLINE  
2025SH0010611	US	US Aorta	23-Jun-2025 13:50	 Final Report (ZZ)	SH	NO IMAGES
2025SH0010289	CR	XR Abdomen - One View	19-Jun-2025 11:47	 Final Report (ZZ)	SH	NO IMAGES
2025SH0010523	US	US Aorta	19-Jun-2025 10:16	 Final Report (ZZ)	SH	ONLINE  
2025SH0010250	CT	CT Left Shoulder and Thoracic Inlet	18-Jun-2025 20:46	 Final Report (ZZ)	SH	ONLINE  

2. If you click on the report or image you want to view - you'll see this preview screen, to view image double-click on the image icon to open up the viewer

Refine your search with Advanced Search

Kestral, Hamish (test)

DOB: 01-Jan-2001 (25 Y) Gender: Male Patient ID: PUA0569


US Right Shoulder	
Priority	Routine
Accession Number	2026AH0003503
Referring Physician	Minnie Mouse
Organization	AH
Status	Final Report (ZZ)
Study Date	12-Aug-2026 10:42
Radiologist	Marshall, Leanne
Location	AH/AH^^^AH^^^^Canopy Imaging - Ahuriri

Final Report (ZZ)

See all patient exams (125)

US
(1 image)

3. You can view previous studies for the patient by choosing studies in the lefthand list. To go back to the previous and click back to the preview screen by hitting the 'x'



- 


referrersupport@canopyimaging.co.nz or DDI: 09 487 2555.

Step 3A – If the break glass symbol is displayed, this means you’re required to follow the below steps to access images and reports

-  Break Glass

I understand that I am about to access information for a patient for whom I am not the primary referring physician on record. I acknowledge that I am doing so only for medical reasons and for the sole benefit of the patient, and in doing so am bound to keep this information private and confidential.

- Please read and acknowledge by ticking the required box before proceeding
- Another search window will appear



Access restricted studies

This feature allows you to view studies for which you were not granted automatic access.

I understand that I am about to access information for a patient for whom I am not the primary referring physician on record. I acknowledge that I am doing so only for medical reasons and for the sole benefit of the patient, and in doing so am bound to keep this information private and confidential.

Note: This feature allows you to view only studies currently available on the PACS. You do not have access to studies for this patient that are acquired at a later date.

Please allow a few minutes for the studies to become available.

☒ I acknowledge and accept the terms of this confidentiality agreement.

Provide the full patient ID

Or


Provide the full patient name and the date of birth

Last name, first name YYYYMMDD

Check Patient

Cancel

- Provide the patient information you have in the relevant fields
- 'Check Patient'
- You will then be shown a patient applicable to your search as below:



A patient that matches your search criteria was found:

HAMISH (TEST) KESTRAL (Male - 25 Y)

Access Patient Studies

Cancel

- Click Access Patient studies
- A list of previous imaging held in Canopy Imaging PACS will appear

If you require remote support to set up IntelConnect, please contact our PACS Team by email: referrersupport@canopyimaging.co.nz or DDI: 09 487 2555.